Immobilization corset

Dear customer, we thank you for the trust placed in an EMO product. By choosing an orthosis from the Corsés Rígidos line you have acquired a quality product and high medical standard. Read these instructions carefully. If you have any questions, please contact your doctor or your specialized establishment.

EFFECTS

Knight corsets: Dislocation. Control of movement in the lumbosacral zone in the antero-posterior, lateral and rotational planes. Taylor-Knight corsets: Dislocation. Control of movement in the lumbosacral and dorsolumbar zones in the antero-posterior, lateral and rotational planes. CLS710 and CLS720, place the posterior valve at the correct height, place the anterior valve facing the anterior one so that both are overlapped on the sides. In CLS710 pass the straps through the buckles on both sides to adjust later. On CLS720 engage the locking wheel BOA, press the wheel until you hear the "click" and turn to adjust, pull the wheel to release. Adjust the cut to the patient's anatomy. It is possible to heat to conform the material locally.

NORMATIVE

The medical device complies with Regulation (EU) 2017/745 of the European Parliament and of the Council on medical devices and the respective national regulations. Class I medical device, non-sterile, without measurement function.

DESCRIPTION AND CHARACTERISTICS

Taylor-Knight Corsets: passive dorsolumbar orthosis formed by an anterior abdominal support attached in its posterior part to a rigid structure, that extends until the first dorsal vertebrae. The models depend on the material of the structure, aluminum CDL508, aluminum pendulum pendulum CDL509, or thermoplastic CDL558. Knight Corset is a lumbosacral passive orthosis formed by an anterior abdominal support attached at the back to a rigid structure. The models depend on the material of the structure, aluminum CLSS01, aluminum pendulum abdomen CLSS05, or thermoplastic CLSS51.

INDICATIONS

Knight and Taylor-Knight corsets:

Vertebral coining. Lumbar trauma.

Hernia. Discal degenerations.

Osteoporosis. Vertebral fractures.

Lumbalgias and dorsalgias that require a temporary immobilization.


Spondylitis. Spondyloarthrosis. After arthroscopic surgery. Postoperative. CLS710 and CLS720: subacute and chronic back pain and back pain, canal stenosis, spondylosis, osteoporosis, degenerative discopathy, degenerative processes, muscle hypotonia, osteoporotic fracture, post-operative rehabilitation, post-operative. CLD718, CLD728 and CLD788, dorsolumbar bivalve corset, bivalve corset made of low density polyethylene. Lateral overlap between 5 cm. The flexibility of the material allows adaptation to the morphology of most patients. The strength of the material maintains the stabilization and immobilization necessary for the pathologies for which it has been designed. CLS710 with strap closures and CLS720 with closure by BOA system. CLD718, CLD728 and CLD788, dorsolumbar bivalve corset, bivalve corset made of low density polyethylene. Lateral overlap between 5 cm. The flexibility of the material allows adaptation to the morphology of most patients. The strength of the material maintains the stabilization and immobilization necessary for the pathologies for which it has been designed. Scapulaxillary straps that facilitate the retroplusion of the shoulders. CLD718 with belt locks and CLD728 with closure by BOA system. CLD788 with belt loops and adjustable height chest support with shoulder straps with adjustable hooks. OVL15 Module for corset 15°. Prefabricated module for the manufacture of stabilization or immobilization corset, with overlap on the front side, made of low density polyethylene. The flexibility of the material allows adaptation to most patients. The strength of the material maintains the stabilization and immobilization necessary for the pathologies for which it has been designed. Possibility of confection with interior lining. Possibility of customization with transfer paper scenery.

In the case of Taylor-Knight, the straps pass over the shoulders and under the arms to the back, where they cross to the front of the waist to be fixed in the corresponding buckles.

CLS710 and CLS720, place the posterior valve at the correct height, place the anterior valve facing the anterior one so that both are overlapped on the sides. In CLS710 pass the straps through the buckles on both sides to adjust later. On CLS720 engage the locking wheel BOA, press the wheel until you hear the "click" and turn to adjust, pull the wheel to release. Adjust the cut to the patient's anatomy. It is possible to heat to conform the material locally.

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QUALITY CERTIFICATES

All the materials used have been tested, comply with the applicable European regulations, are of high quality and offer maximum safety and comfort during use. Oko-Tex Standard 100 Certified: Textile eco-label awarded following the completion of a catalog of test criteria ensuring the absence of substances harmful to health.

Certifies the absence of harmful substances in all materials with which a garment is made.

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CLS710 et CLS720, lumbosacral bivalve corset: bivalve corset made of low density polyethylene, with lateral overlap between 5 cm. The flexibility of the material allows adaptation to the morphology of most patients. The strength of the material maintains the stabilization and immobilization necessary for the pathologies for which it has been designed. CLS710 with strap closures and CLS720 with closure by BOA system. CLD718, CLD728 and CLD788, dorsolumbar bivalve corset, bivalve corset made of low density polyethylene. Lateral overlap between 5 cm. The flexibility of the material allows adaptation to the morphology of most patients. The strength of the material maintains the stabilization and immobilization necessary for the pathologies for which it has been designed. Scapulaxillary straps that facilitate the retroplusion of the shoulders. CLD718 with belt locks and CLD728 with belt loops and adjustable height chest support with shoulder straps with adjustable hooks. OVL15 Module for corset 15°. Prefabricated module for the manufacture of stabilization or immobilization corset, with overlap on the front side, made of low density polyethylene. The flexibility of the material allows adaptation to most patients. The strength of the material maintains the stabilization and immobilization necessary for the pathologies for which it has been designed. Possibility of confection with interior lining. Possibility of customization with transfer paper scenery.

Textile materials: Hand wash with warm water (maximum 30 °) and neutral soap. Dry at room temperature. Do not store or iron the garment. Contact with ointments, creams, alcohols and solvents may deteriorate the garment. Plastic and metal parts: Wash with a cloth moistened with water, do not use solvents or other chemicals. Dry mainly areas containing metal parts.

CLD718, CLD728 and CLD788, place the posterior valve at the correct height, place the anterior valve facing the anterior one so that both are overlapped on the sides. In CLD718 pass the straps through the buckles on both sides to adjust later. On CLD728 engage the locking wheel BOA, press the wheel until you hear the "click" and turn to adjust, pull the wheel to release. Adjust the cut to the patient's anatomy. It is possible to heat to conform the material locally.

INSTRUCTIONS FOR USE AND PLACEMENT

Knight and Taylor-Knight corsets: Place the corset wrapping the trunk, adapting it to the patient's pelvic contour with the adjustment straps fastened at the front. The frame should support your lower buttocks, as low as possible, preventing you from bothering when sitting on a chair, and on the upper back. Modify the structure if necessary to adapt it to the patient.

WASHING AND MAINTENANCE

Wash with a damp cloth and mild soap. During washing, if the garment has Velcro, it must be glued together or kept closed. Dry at room temperature. Do not expose the brace directly to a direct heat source. Store at room temperature.

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