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IF020 Splint Stack In Polypropylene
IF030 Splint Stack In Polyethylene
The Stack splint is designed to achieve the immobilization of the distal interphalangeal joint in every finger. This orthosis is indicated for painful stages of the joint, usually caused by direct trauma.
PLACEMENT: Enter your fingertip affected by the brace to reach its bottom, the maximum coupling to the finger brace. Surround your finger with the strap attached or riveted to the splint in proximal part, adjusting without exerting excessive pressure. Verify that the bending of the proximal interphalangeal joint of the finger is not restricted. If this has been limited, see your prosthetist will proceed to cut the brace properly. In the event that the brace is not perfectly coupled to the morphology of the finger, the prosthetist can proceed to manual shaping or by applying hot air to the IF030 (polyethylene) or by immersing the splint in hot water for IF020 (polypropylene.)
IF020 is supplied in bags of 10 pcs. without reference IF20B tapes.
SIZES: 1, 2, 3, 4, 5, 6, 7, 8

Restraint In Thermoplastic Interphalangeal IF010
Indicated in cases of immobilization and rest of the interphalangeal joints, painful stages usually produced as a result of direct trauma.
PLACEMENT: Cover the palm area of the affected finger without any gaps, forming manually if necessary, attach the velcro on the thermoplastic adhesive at the joints to close the splint with rubber bands.
SIZE (LENGTH OF THE FINGER MEASURED BY THE VOLAR IN CM)
1 (6.3), 2 (6.8), 3 (7.1), 4 (7.5), 5 (8.5), 6 (9.5), 7 (10), 8 (10.3), 9 (10.8).

Immobilizer IF040 IFD "Rana Federica"
Indicated in cases of immobilization and rest of the interphalangeal joints, painful stages usually produced as a result of direct trauma.
PLACEMENT: Cover the palm area of the affected finger, closing with a double curve, the side flaps over the middle and proximal phalanges. The splint can be cut

with shears, then going over the edges to avoid sharp edges. Repetitive bending of the aluminum produces the rupture.
SIZE (MAXIMUM LENGTH OF THE FINGER IN CM.) S (6.3), M (7.1) L (7.7)
IF050 Ferula Progressive Finger Extension
Shown in Dupuytren's contracture, a disease of chronic and progressive course which leads to shortening and thickening of the fascial strips on the palmar fascia, with consequent bending deformation in proximal interphalangeal joint, metacarpophalangeal joint, or both. The application of this splint gets slow deformation that causes this disease, acting as a progressive correction of flexion of the fingers.
PLACEMENT: Hosting the affected finger palmar supporting the workpiece rigid plastic surround with tape so that the foam pad is placed on the dorsal proximal interphalangeal joint, adjust the belt with the buckle to length to hold the finger allowing progression through the crossbow bolt, allowing minimal and progressive increases of the pressure required for correction. To correctly position the splint must match the convexity of the splint with the angle between the proximal phalanx and palm, leaving the metal rivet closer to the palm.
SIZES: S, M, L

Flash and Precautions
Adjust with the tapes does not require excessive pressure.
To avoid unwanted rigidities, the splint should not be used longer than directed by the physician.
Do not expose the brace to a flame or heat source at high temperatures, it may deform or damage it.
Washable with mild soap and water. Air dry, or towel.

Regulations
Regulations The medical device complies with Regulation (EU) 2017/745 of the European Parliament and of the Council on medical devices and the respective national regulations. Class I medical device, non-sterile, without measurement function.

Precautions
Follow the manufacturer's directions for proper orthosis placement. Follow the instructions or recommendations of the prescriber or healthcare personnel. For a personalized adaptation, consult a specialist. If there is ongoing pain, it is advisable to contact a healthcare professional.
If an allergic reaction should occur remove the product and consult with the prescribing physician. In case of minor annoyances caused by excessive sweating we recommend the use of a tubular cotton. Positioned such that compression is not excessive. Do not expose the garment to a flame or heat source at very high temperatures. Store at room temperature, preferably in the original packaging. Guarantee within the framework of laws.

Notice
To dispose of the product, use the solid waste container. For the disposal of the container, depending on its material, use the paper / cardboard container or the plastic container. In any case, you must strictly comply with the regulations of the country. Any serious incident related to the product must be reported to the manufacturer and to the competent authority of the Member State in which the user and / or the patient are established.

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